

THREE IN ONE *Entonox Demand valve Kit*

Pressure Regulator

High-pressure Hose

Demand Valve

Silicone Mask

Features

- Durable anodized aluminum body.
- Demand Valve provides a flow rate of 40LPM and a Manual flow rate up to 160LPM.
- Allows both hands to remain on the mask when depressing the Manual Ventilation Button.

It's to be used for patients (adult and child) greater than 40kg who need first aid, and/or oxygen and N2O therapy.

Features




Oxygen supply :

- Inlet Connection: International standard yoke (Inlet pressure around 2900 psi, without flow adjust knob)
- Output Pressure: 60 PSI
- Flow rate: 0-25 LPM



Demand valve :

- Gas inlet connection : 5/8" DISS
- Gas outlet : 22mm*15mm
- Pressure : 40~60PSIG
- Flow rate : Provides 40LPM Manual Flow Rate and up to 160LPM flow rate on Demand by the patient.
- Pressure Relief Valve: 60cm H2O Features



FDA ISO 13485:2003 GMP CE 0434