



OXYLITRE “petite elite” PORTABLE SUCTION PUMP HIGH SUCTION

Specification

- Low Vacuum High Flow
- Single 1000ml Disposable SERRES Jar & Liner System
- Robust High Impact Resistant Moulded Enclosure
- Encapsulated Bacterial/Hydrophobic Filter Disk
- Luminous Touch Panel Control
- Instant Maximum Vacuum Override Button
- Battery Level/ Charging Indicator
- 10 to 150mmHg.
- Flowrate >25Lpm
- 110volt – 240volt AC
- 12volt – 28volt DC
- 12 volt DC Rechargeable Battery, 4.3mAh
- Run time >1hr 20min
- Weight 4.9Kg
- Approximate Size 310mm x 305mm x 170mm
- Peak Sound Output at MAX Vacuum from 1m: 65 decibels